

WO 2004/003520

2/PRTS.

10/519340

PCT/GB2003/002775  
DT09 Rec'd PCT/PTO 22 DEC 2004

- 1 -

## TWO-STAGE PARTICLE-SIZE ANALYSER

5 The present invention relates to a method and apparatus for collecting particles suspended in a fluid.

10 Collecting the particulate matter suspended in a fluid, e.g. air, is an important stage of air quality assessment, atmospheric science and aerosol technology, and particles collected from a fluid are analysed by various chemical and physical methods for particulate matter characterisation.

15 There are two methods of particle characterisation (i) bulk analysis and (ii) size selective analysis. The latter usually involves describing the particle size distributions and so the size selective collection of particles is an important stage in their characterisation and the present invention relates to the size selection of particles.

20 A known method of characterising aerosol particles size distributions is based on the deposition of particles onto substrates in a cascade impactor and further analysis of the deposits (e.g. by gravimetric or chemical analysis). In a cascade impactor particles of different sizes are collected onto different substrates due to the difference in their inertia. The selectivity of deposition is achieved by means of a number of air jets with specific aerodynamic characteristics. Each stage of an impactor has a different jet facing the substrate where particles are collected. Thus, an impactor enables a set of mass concentrations in various size ranges (size sections) to be  
25 obtained.

Usually a cascade impactor is used to collect particles larger than 0.3  $\mu\text{m}$ . This method has a limitation and it is very difficult to apply it to particles smaller than 0.3  $\mu\text{m}$ .

30

WO 2004/003520

PCT/GB2003/002775

- 2 -

Another method relevant to obtaining aerosol size distributions is the deposition of aerosol particles onto a fibre or membrane filter. In this method a size selective inlet is often used to remove particles larger than a certain size, e.g. 10  $\mu\text{m}$ . All particles passed through the inlet are collected onto a filter and are analysed later. This method is simpler to use than cascade impactors. Various size selective inlets are used along with a filter to sample the mass fraction of an aerosol; for instance  $\text{PM}_{10}$ ,  $\text{PM}_{2.5}$  or  $\text{PM}_1$  (where the figure indicates the cut off aerodynamic diameter of the inlet). The filter method enables a wide range of particles to be collected, even particles smaller than 0.3  $\mu\text{m}$ . However it has a limited capability for obtaining information about particle sizes and, in particular, the major drawback of this method is its inability to deliver the size resolved information so size distributions of the particulate matter cannot be obtained with this technique.

We have devised an improved method and apparatus for collecting particles from a fluid.

According to the invention there is provided a method for selective deposition of suspended particles from a fluid which method comprises passing the fluid sequentially over a first collector adapted to collect larger particles and then over a second collector adapted to collect smaller particles, which second collector comprises a chamber in which there is at least one net or another material containing fibres placed across the chamber.

The invention also provides a particle collector for collecting and sampling particles in a fluid which comprises sequentially (i) an inlet, (ii) a first collector adapted to collect larger particles and (iii) a second collector adapted to collect smaller particles comprising a chamber in which there is at least one net placed across the chamber and a flows means able to sustain a flow of fluid sequentially through the inlet, first collector and second collector.

30

WO 2004/003520

PCT/GB2003/002775

- 3 -

The nets can be any structure which has the equivalent effect to nets, e.g. can be woven, knitted or formed of fibres so that the effect is similar to nets in removing particles; for example they can also be rigid or semi rigid.

- 5 By larger particles is meant particles larger than those collected in the second collector. In general this will mean particles of sizes above about 0.3  $\mu\text{m}$ .

There optionally can be further collecting media through which the fluid subsequently will flow in use.

10

Preferably the first collector comprises a cascade impactor or a sedimentation cell, e.g. containing set of parallel horizontal partitions.

- 15 The second collection collector can comprise at least one net and preferably at least two nets of different mesh sizes mounted within a container, so that the fluid passes sequentially through the nets. There can be three, four, five or more nets.

- 20 Thus, particles of different sizes are collected on different nets. In a four net construction for example the first net faces the flow and collects the largest particles, e.g. greater than 100 nm; the particles smaller than 10 nm penetrate through the first net; the second net collects the particles in the size range from 10 to 30 nm; the particles smaller than 30 nm but larger than 10nm penetrate through the second net and they are collected by the third net; particles smaller than 10nm (e.g. from 1 to 10 nm) are collected by the fourth net.

25

The nets can be either identical or different. Different nets can be used to increase the size range of particles to be deposited. For example the first net can have a mesh opening of 120  $\mu\text{m}$ ; the second net can have a mesh opening of 40  $\mu\text{m}$ ; the third net can have a mesh opening of 20  $\mu\text{m}$  and the fourth net can have a mesh opening of 10  $\mu\text{m}$ .

30

WO 2004/003520

PCT/GB2003/002775

- 4 -

The net sampling is applicable only for sub-micron size particles, for instance for particles smaller than about 0.3  $\mu\text{m}$ . So the first collecting collector preferably collects particles above this size.

5

When a sedimentation cell with a plurality of partitions is used as the first collector, the particles are separated due to gravitational settling onto the partitions and can be analysed later. The sedimentation of particles depends on their size. Thus analysing the different parts of the pile of partitions make it possible to obtain additional information about the size distribution of the particles.

10

Preferably the collection of particles on the net(s) takes place at controlled humidity and preferably there is a humidity control unit which is incorporated between the inlet and the large particle collector.

15

The invention is suitable for use with aerosols and, in use with an aerosol, the aerosol particles are introduced into the inlet and after that go into the first section of the first collector (e.g. the first stage of the cascade impactor). A fraction of particles of the higher collection ability is collected by the first stage. The rest of the particles goes further with the flow and is deposited onto the next stages. Every stage collects particles of certain sizes. After passing all the stages of the cascade impactor, the flow goes into the net sampler where smaller particles are deposited according to their efficiency. A fraction of particles of the higher collection ability is collected by the first net. The rest of particles goes further and is collected by the next nets.

20

25

The invention enables there to be delivered size resolved information so a much more accurate method size distribution of the particulate matter can be obtained with this technique and it was very surprising that the combination of the two different separators gives such improved results.

30

WO 2004/003520

PCT/GB2003/002775

- 5 -

The invention is illustrated in the accompanying drawings in which

Fig. 1 shows an existing collector,

Fig. 2 shows schematically a net collector which can be used,

- 5 Fig. 3 shows schematically a simple collector according to the invention and  
Fig. 4 shows schematically a more detailed collector according to the invention.

Referring to fig. 1, existing collectors for use with aerosols comprise a size selective  
preseparator (21) (e.g. a cyclone with 10  $\mu\text{m}$  cut off aerodynamic size), inlet (22),  
10 filter (23) and outlet (24).

In use the aerosol passes through cyclone preseparator (21) which collects particles  
above 10  $\mu\text{m}$  and the aerosol particles pass through inlet (22) and are then deposited  
onto a fibre or membrane filter (23) and the air then passes out through outlet (24).  
15 All the particles which have passed through the inlet are collected onto the filter and  
are analysed later.

Referring to fig. 2 which shows a net collector, there is a container (16) with inlet  
(11), nets (12) and outlet (13). In use particles of different sizes are collected on  
20 different nets. Seven nets are shown for illustration with the first net facing the flow,  
but in a four net construction; for example the first net faces the flow and collects the  
largest particles e.g. greater than 100 nm; the particles smaller than 10 nm penetrate  
through the first net; the second net collects the particles in the size range from 10 to  
30 nm; the particles smaller than 30 nm but larger than 10nm penetrate through the  
25 second net and they are collected by the third net; particles smaller than 10nm (e.g.  
from 1 to 10 nm) are collected by the fourth net. The particles can be removed from  
the net and analysed which makes it possible to obtain additional information about  
the size distribution of aerosol particles.

WO 2004/003520

PCT/GB2003/002775

- 6 -

Referring to fig. 3 there is a sedimentation cell (17) containing horizontal partitions (15) and inlets and outlets (14), (11), and (13).

5 In use the particles are separated in separator (17) due to gravitational settling onto the partitions (15) and these particles can be analysed later. The sedimentation of particles depends on their size. Thus analysing the different parts of the pile of partitions makes it possible to obtain additional information about the size distribution of aerosol particles. After leaving (17) the fluid passes through inlet (11) to net separator (16) which functions are as described in fig. 2.

10

Referring to fig. 4, this illustrates a wide range of aerosol samplers designed to collect selectively aerosol particles in a wide range of sizes from 1nm to 30  $\mu\text{m}$  aerodynamic diameter under a constant controlled humidity. The flow rate is from 1 to 30 l/min and the sampling humidity (inside the sampling system) is from 30 to 95%.

15

The collector and sampling system consists of a net sampler (1), humidity control unit (2), cascade impactor (3), aerosol chamber (4), inlet (5), flow meter (6), saturator (7), pump (8) and outlet (9) with aerosol filter.

20 In use this system provides sampling under a constant humidity that could be set using the humidity control unit button on humidity control unit (2). Aerosol enters the saturator (7) through the inlet (5). After the saturator the aerosol goes into the aerosol chamber (4) where water vapour condenses onto particles. The chamber (4) is connected to the humidity control unit (2). If humidity is lower than required the  
25 heater in the saturator is turned on by the humidity controller. It gives more water vapour and humidity is increased.

After the humidity control unit (2) aerosol enters the cascade impactor (3) where particles larger than 0.25  $\mu\text{m}$  are collected onto impactor slides: Microscope Slides  
30 (Size 76 x 26 mm; thickness 1.0 -1.2 mm). The particles smaller than 0.25  $\mu\text{m}$  are

WO 2004/003520

PCT/GB2003/002775

- 7 -

then collected by the net sampler (1). The particles can then be analysed according to their size and the size distribution is calculated.

The cascade impactor

5

The size bands of a cascade impactor are influenced by the flow rate. At the flow rate 20 l/min 50% particle retention efficiency, aerodynamic diameters are shown in Table 1.

10 Table 1. 50% Particle Retention Efficiency aerodynamic diameters for May cascade impactor used in the prototype

Impactor stage number	1	2	3	4	5	6	7
50% Particle retention efficiency aerodynamic diameter, $\mu\text{m}$	20	8	4	2	1	0.5	0.25

15 The net sampler

There are two options shown as examples:

- (i) the basic configuration of the net sampler with 4 size sections and
- (ii) 8-section net sampler.

20

At the flow rate 20 l/min maximal collection efficiency aerodynamic diameters for the basic configuration of the net sampler are shown in Table 2.

25

WO 2004/003520

PCT/GB2003/002775

- 8 -

Table 2. Maximal collection efficiency aerodynamic diameters for the basic configuration of the net sampler

Net sampler section number	1	2	3	4
Maximal collection efficiency aerodynamic diameter, nm and ( $\mu\text{m}$ )	128 (0.128)	32 (0.032)	8 (0.008)	2 (0.002)

5

Maximal collection efficiency aerodynamic diameters (at the flow rate 20 l/min) for 8-section net sampler are shown in Table 3.

Table 3. Maximal collection efficiency aerodynamic diameters for 8-section net sampler

10

Net sampler section number	1	2	3	4	5	6	7	8
Maximal collection efficiency aerodynamic diameter, nm	128	64	32	16	8	4	2	1

The deposits on nets can be analysed separately. The size distribution of an aerosol is determined from chemical analysis or gravimetric measurements.

15